



Martin-Luther-Universität Halle-Wittenberg
 Institut für Geowissenschaften
 FG Mineralogie/Geochemie



Prof. Dr. Dr. H. Pöllmann
 Von-Seckendorff-Platz 3, 06120 Halle,
 Tel.: +49/345/5526110 Fax: +49/0345/5527180
 e-mail: herbert.poellmann@geo.uni-halle.de
 Ort: H3 0.30.1

PANalytical EMPYREAN


Computed Tomography (CT)

The PANalytical EMPYREAN X-ray diffractometer is equipped with a Mo-tube (/Cu-tube), a rotational sample stage for CT and a state of the art detection system, the GaliPIX^{3D}, that is capable of collecting frames suitable for 3D processing in computed tomography.

CT offers the possibility to **analyze and visualize different features** (e.g. amount/type/connectivity of pores, cracks, inclusions, mineral distribution, preferred orientations, wall-thickness, surface area) **in a fast and non-destructive way through 3-dimensional virtual reconstructions of the sample.**

Applications: Building materials, ceramics, composite compounds, foams, electronics, rocks.



Experimental Setup - Overview



- 1) Mo- or Cu-tube in point focus
- 2) incident beam optic & filter
- 3) CT rotational sample stage
- 4) sample holder with sample
- 5) GaliPIX^{3D} detector

Experimental parameters (anode material, voltage/current, filter materials, rotation range, step size, time per step, PHD levels) can be adjusted according to the sample.

Example: Autoclaved Aerated Concrete



The **3D pore-size(-volume) distribution** (which has an influence on strength, permeability, density, thermal insulation, etc.) can be determined and visualized with CT on a laboratory X-ray diffractometer.